



MK1 Attachment Dental Prosthesis Silver Partial Denture Precision Attachments

Our Product Introduction

Basic Information

- Place of Origin: Shenzhen, China
- Brand Name: ADL
- Certification: ISO/CE/FDA
- Minimum Order Quantity: Negotiable
- Price: Negotiable
- Packaging Details: Plastic bags and box
- Delivery Time: 72 Hours
- Payment Terms: Paypal or T/T
- Supply Ability: 630 Technicians



Product Specification

- Durability: Highly Durable
- Precision: High Precision
- Usage: Dental Prosthesis
- Color: Silver
- Installation: Easy To Install
- Adjustability: Adjustable For Different Prosthesis Sizes
- Product Type: Precision Attachment
- Function: Precision Attachment For Dental Prosthesis
- Highlight: **MK1 Attachment Dental Prosthesis Silver ,
MK1 Dental Prosthesis Silver ,
MK1 Dental Prosthesis Silver Partial Denture**



for more products please visit us on adidentalab.com

Product Description

Precision Silver Easy To Install MK1 Dental Prosthesis Attachment

Product Description:

The MK1 Precision Attachment is a product that is highly durable, with the ability to withstand the wear and tear of regular use. This product is designed to be long-lasting, ensuring that it can be used for a long time without the need for frequent replacement. The MK1 Precision Attachment is a precision attachment for dental prostheses that is designed to be used in a variety of applications. This product is ideal for use in both fixed and removable dental prostheses, providing high-precision results that are both accurate and consistent.

The MK1 Precision Attachment is a precision attachment that is designed to be highly precise, ensuring that the results obtained are accurate and consistent. This product is designed to be easy to use, ensuring that even those who are new to dental prosthetic applications can use it with ease.

The MK1 Precision Attachment is a highly precise product that is ideal for use in a variety of dental prosthetic applications. This product is designed to be both durable and long-lasting, ensuring that it can be used for a long time without the need for frequent replacement.

The MK1 Dental Attachment is a precision attachment that is designed to be used in a variety of dental prosthetic applications. This product is highly precise, ensuring that the results obtained are accurate and consistent. The MK1 Dental Attachment is also highly durable, ensuring that it can be used for a long time without the need for frequent replacement.

In conclusion, the MK1 Precision Attachment Denture is a highly precise product that is ideal for use in a variety of dental prosthetic applications. This precision attachment is designed to be both durable and long-lasting, ensuring that it can be used for a long time without the need for frequent replacement. The MK1 Dental Attachment is highly precise and accurate, providing consistent results that are both accurate and reliable.

Features:

Product Name: MK1 Precision Attachment

Installation: Easy To Install

Precision: High Precision

Package Contents:

1 MK1 Precision Attachment

1 Installation Tool

Adjustability: Adjustable For Different Prosthesis Sizes

Function: Precision Attachment For Dental Prosthesis

Technical Parameters:

Adjustability	Adjustable For Different Prosthesis Sizes
Compatibility	Universal
Function	Precision Attachment For Dental Prosthesis
Color	Silver
Durability	Highly Durable
Installation	Easy To Install
Usage	Dental Prosthesis
Precision	High Precision
Package Contents	1 MK1 Precision Attachment, 1 Installation Tool
Product Type	Precision Attachment

Applications:

Introducing the MK1 Precision Attachment, a universal and durable precision attachment for dental prosthesis manufactured by ADL in Shenzhen, China. This product has been certified by ISO, CE, and FDA, ensuring its high quality and safety for dental use.

The MK1 Precision Attachment is compatible with various dental prosthesis and is easy to install with the included installation tool. It is highly durable, ensuring longevity and reliability for its users.

The MK1 Precision Attachment is ideal for patients who require partial dentures and need a reliable and precise attachment. It is suitable for various occasions and scenarios, such as:

Patients who need partial dentures due to tooth loss or damage

Dentists who need a reliable and durable precision attachment for their patients' dental prosthesis

Dental labs who require a universal precision attachment for their dental prosthesis production

The MK1 Precision Attachment is available for purchase with negotiable price and minimum order quantity. It is packaged in plastic bags and box for safe and easy handling and delivery time is within 72 hours. Payment options include Paypal or T/T, and ADL has a supply ability of 630 technicians to ensure timely delivery.

With the MK1 Precision Attachment, patients can have a comfortable and secure fit for their partial dentures, while dentists and dental labs can ensure a precise and reliable attachment for their patients' dental prosthesis. Get your MK1 Precision Attachment today for a

better dental experience.



Advanced Dental Laboratory



+86 18665847934



info@adidentalab.com



adidentalab.com

Plant 1, Fuqiao District 2, Qiaotou Community, Fuhai Street, Baoan District, Shenzhen City, Guangdong Province, P.R.China